

VC211 EXPERIMENT E5 DATASHEET: PRECIPITATION & WATER PURITY										
STDNT:		ID:		SECTION#:		TA:				
GRP#:										
PART 1: What is a precipitate. Each 2 students tests 1 sample once										
	CuSO ₄ 5mL (0.1M) Color	BaCl ₂ 5mL (0.1M) Color	CuSO ₄ + BaCl ₂ Color	ppt↓ Yes/No	Filtrate observed property					
2students										
2students										
PART 2.A. Is Precipitation Predictable? Group efforts, test only Group I or Group II as assigned, 2 raw reactions/student. Clear = means no precipitate, then record solution color.										
CATION GROUP I	CATION GROUP II	REF WATER	Cl ⁻ 2drops	CrO ₄ ²⁻ 2drops	I ⁻ 2drops	C ₂ O ₄ ²⁻ 2drops	S ²⁻ 2drops	SO ₄ ²⁻ 2drops	SPECTATOR IONS	
Cations no. drops→		2drops	2drops	2drops	2drops	2drops	2drops	2drops	GROUP I	GROUP II
Na ⁺	K ⁺	clear / colorless								
Ba ²⁺	Mn ²⁺	clear / colorless								
Mg ²⁺	Ca ²⁺	clear / colorless								
Co ²⁺	Sr ²⁺	clear / colorless								
Ni ²⁺	Cr ³⁺	clear / colorless								
Cu ²⁺	Fe ³⁺	clear / colorless								
Al ³⁺	Zn ²⁺	clear / colorless								
Pb ²⁺	Ag ⁺	clear / colorless	White ppt↓	Brown ppt↓	Yellow ppt↓	White ppt↓	Black ppt↓	White ppt↓		
PART 3. Conce. & Precip.: Each team uses Table 3 & design different reactions than Table, minimum 2 reactions per student										
TABLE 5: RECORD YOUR RAW DATA HERE similar to Table 4. Add only 2 drops of each reactant				TABLE 3: DESIGN REACTIONS FROM HERE			TABLE 4: SAMPLE REACTIONS DESIGN			
REACTION #	REACTANT #1 & CONC	REACTANT #2 & CONC	OBSERVATIONS	REACTION #	REACTANT #1	REACTANT #2	Reactant #	REACTANT #1 & Conc.	REACTANT #2 & Conc.	OBSERV.
				I	Pb(NO ₃) ₂	KI	II-1	0.10M Pb(NO ₃) ₂	0.10M NaOH	
				II	Pb(NO ₃) ₂	NaOH	II-2	0.10M Pb(NO ₃) ₂	1.0M NaOH	
				III	AgNO ₃	KI	II-3	0.01M Pb(NO ₃) ₂	0.01M NaOH	
				IV	ZnSO ₄	NaOH	IV-1	0.10M ZnSO ₄	0.10M NaOH	
				V	CaCl ₂	K ₂ C ₂ O ₄	IV-2	0.10M ZnSO ₄	1.0M NaOH	
				VI	CaCl ₂	NaOH	IV-3	0.01M ZnSO ₄	0.01M NaOH	
							V-1	0.10M CaCl ₂	0.10M K ₂ C ₂ O ₄	
							V-2	0.01M CaCl ₂	0.01M K ₂ C ₂ O ₄	
PART 4. Solvent Pollution & Precip.: Group efforts, each 2 students study solubility of 1 solid										
TABLE 6: Solids Solubility in Polar & Non-Polar Solvents (total 6 samples to test)										
Solid Type (thoroughly dry inside test tubes)	Ionized Water 2mL	Acetone 2mL	Hexane 2mL	Note: Use the solubility table from CH 4, VC210 to predict if precipitate is formed from mixing the supernatants						
I. CaCl ₂ <0.2g										
II. K ₂ C ₂ O ₄ <0.2g										
Supernatant (I + II)										